


HPLC REPORT

Sample: C-16P
Analyst: WGJ
Lot. No.: P140926-HX608706
Column: SinoChrom ODS-BP, 4.6*250mm, 5µm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
Gradient:
0.0min 16% 84%
25.0min 41% 59%
25.1min 100% 0%
30.0min Stop
Volume: 20µl
Wavelength: 220nm
Flow rate: 1.0ml/min

Analyzed date: 2014-10-10



Rank	Time	Area	Conc.
1	16.537	38219.027	0.9356
2	17.720	141897.844	3.4735
3	18.017	3477532.000	85.1261
4	18.017	99114.883	2.4262
5	21.153	304926.094	7.4643
6	25.317	23462.865	0.5743

Total 100