


MASS SPECTROMETRY REPORT



Sample Description
 Analyzed date: 2014/10/10
 Analyst: Shen
 Sample: C-16P
 M.W.: 1995.32
 Lot. No.: P140926-HX608708

Instrument	SHIMADZU LCMS-2020		
Probe:	ESI	Probe Bias:	+ 4.5kv
Nebulizer Gas Flow:	1.5L/min	Detector:	1.2kv
CDL:	-20.0v	T. Flow:	0.2ml/min
CDL Temp.:	250 °C	B. Conc.:	50% H ₂ O/50% ACN
Block Temp.:	400 °C		