

MASS SPECTROMETRY REPORT


30-Sep-2017
09:23:45

Probe: ESI
Cone: 50v
Desolvation Temp :350

Capillary:3.00KV
Extractor: 5v
Gas Flow : 350

Scan ES+
2.28e5

170930-9 9 (0.166)



Sample Description

Analyzed date: 2014-09-30
Analyst: YU
Sample: C-16
M.W.: 1351.41
Lot. No.: P140921-HS608567

Instrument

Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp.: 250 °C
Block Temp.: 200 °C

Waters ZQ2000

Probe Bias: + 4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. Conc.: 50%H2O/50%ACN