

HPLC REPORT

Sample: K-12

Analyzed date: 2014-10-17

Analyst: ZK

Lot. No.: P141010-WD608789

Column: Kromasil 100-5C18, 4.6*250mm, 5µm

Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile

Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water

Gradient :

	A	B
0.0min	18%	82%
25.0min	43%	57%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl

Wavelength: 220nm

Flow rate: 1.0ml/min

mV

-220

-200

-180

-160

-140

-120

-100

-80

-60

-40

-20

0

-20

2

4

6

8

10

12

14

16

18

20

22

24

26

min

Rank	Time	Quantity	Area	Height
1	10.293	6.9111	204721	15914
2	10.601	90.5552	2682436	173891
3	11.202	1.9235	56979	5740
4	12.096	0.6102	18075	2090
Total		100	2962211	197635