


HPLC REPORT

Sample: T-14
Analyst: WGJ
Lot. No.: P160922-YJ608754
Column: SinoChrom ODS-BP, 4.6*250mm, 5µm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
Gradient:
0.0min 16% 84%
25.0min 41% 59%
25.1min 100% 0%
30.0min Stop
Volume: 20µl
Wavelength: 220nm
Flow rate: 1.0ml/min

Analyzed date: 2016-10-25



Rank	Time	Area	Conc.
1	10.788	173990.297	5.7699
2	10.995	2588881.250	85.8528
3	11.297	252617.422	8.3773
Total		100	