

HPLC REPORT

Sample: C-13
Analyst: WGJ
Lot. No.: P140922-YJ608766
Column: SinoChrom ODS-BP, 4.6*250mm, 5µm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water

Analyzed date: 2014-09-25

Gradient:


	A	B
0.0min	21%	79%
25.0min	46%	54%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl

Wavelength: 220nm

Flow rate: 1.0ml/min

mV



Rank	Time	Area	Conc.
1	7.653	12927	0.2471
2	8.333	99638	1.904
3	8.476	4889690	93.46
4	8.721	97447	1.863
5	12.828	132179	2.526
Total			100