


HPLC REPORT

Sample: C-14
Analyst: WGJ
Lot. No.: P140922-YJ608772
Column: SinoChrom ODS-BP, 4.6*250mm, 5µm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
Gradient:

	A	B
0.0min	17%	83%
25.0min	42%	58%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl
Wavelength: 220nm
Flow rate: 1.0ml/min

Analyzed date: 2014-09-25



Rank	Time	Area	Conc.
1	10.483	31678.395	1.2621
2	12.782	94384.938	3.7604
3	12.995	2347334.500	93.5194
4	12.995	36599.051	1.4581
Total		100	