


# MASS SPECTROMETRY REPORT



Sample Description  
 Analyzed date: 2014/10/9  
 Analyst: Shen  
 Sample: C-13  
 M.W.: 1529.62  
 Lot. No.: P140922-YJ608782

Instrument	SHIMADZU LCMS-2020		
Probe:	ESI	Probe Bias:	+ 4.5kv
Nebulizer Gas Flow:	1.5L/min	Detector:	1.2kv
CDL:	-20.0v	T. Flow:	0.2ml/min
CDL Temp.:	250 °C	B. Conc.:	50% H <sub>2</sub> O/50% ACN
Block Temp.:	400 °C		