



Sample Description
 Analyzed date: 2014/9/29
 Analyst: Shen
 Sample: C-16
 M.W.: 1915.34
 Lot. No.: P140921-HX608708

Instrument SHIMADZU LCMS-2020
 Probe: ESI Probe Bias: + 4.5kv
 Nebulizer Gas Flow: 1.5L/min Detector: 1.2kv
 CDL: -20.0v T. Flow: 0.2ml/min
 CDL Temp.: 250 °C B. Conc.: 50% H2O/50% ACN
 Block Temp.: 400 °C