

HPLC REPORT


Sample: K-12
 Analyst: WGJ
 Lot. No.: P140705-YJ594475
 Column: SinoChrom ODS-BP, 4.6*250mm, 5µm
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water

Analyzed date: 2014-07-09

Gradient:

	A	B
0.0min	9%	91%
25.0min	34%	66%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl
 Wavelength: 220nm
 Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	10.637	41269	0.92
2	11.307	182228	4.05
3	11.475	4152814	92.23
4	11.809	1259052	2.80
Total			100