


# HPLC REPORT

Sample: T-10P Analyzed date: 2014-7-17  
Analyst: ZK  
Lot. No.: P140710-SY594216  
Column: Kromasil 100-5C18, 4.6\*250mm, 5µm  
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile  
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water  
Gradient :

	A	B
0.0min	15%	85%
25.0min	40%	60%
25.1min	100%	0%
30.0min		Stop

Volume: 20µl  
Wavelength: 220nm  
Flow rate: 1.0ml/min



Rank	Time	Quantity	Area	Height
1	8.438	1.3784	21799	2258
2	9.369	1.7676	27954	5224
3	9.500	96.2020	152143	188506
4	10.198	0.6520	10312	1305
Total		100	1581496	197293