

HPLC REPORT

Sample: C-10 Analyzed date: 2014-06-11

Analyst: YHY

Lot. No.: P140601-YJ587979

Column: SinoChrom ODS-BP, 4.6*250mm, 5µm

Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile

Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water


Gradient:	A	B
0.0min	20%	80%
25.0min	45%	55%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl

Wavelength: 220nm

Flow rate: 1.0ml/min

mV



Rank	Time	Area	Conc.
1	9.613	100347	1.604
2	10.054	5870870	93.83
3	10.267	234276	3.744
4	10.438	51693	0.8261
Total		6257186	100