

# HPLC REPORT

Sample: F-11 Analyzed date: 2014-06-12  
 Analyst: YHY  
 Lot. No.: P140607-YS587985  
 Column: SinoChrom ODS-BP, 4.6\*250mm, 5µm  
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile  
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water  
 Gradient:

	A	B
0.0min	22%	78%
25.0min	47%	53%
25.1min	100%	0%
30.0min		Stop

Volume: 20µl  
 Wavelength: 220nm  
 Flow rate: 1.0ml/min  
 mV



Rank	Time	Area	Conc.
1	15.324	38833	0.1844
2	15.649	64092	0.3043
3	16.226	573448	2.723
4	16.533	492790	2.34
5	17.357	171883	0.8162
6	17.602	18790360	89.22
7	17.938	928833	4.41

Total 21060239 100