


HPLC REPORT

Sample: K-11(P) Analyzed date: 2014-07-16
Analyst: LLR
Lot. No.: P140710-MX594717
Column: Kromasil C18, 4.6*250mm, 5µm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
Gradient:
0.0min 15% 85%
25.0min 40% 60%
25.1min 100% 0%
30.0min Stop
Volume: 20µl
Wavelength: 220nm
Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	11.043	6993	0.334
2	12.356	1894315	90.47
3	12.489	113620	5.426
4	13.126	17778	0.8491
5	14.479	61151	2.921
Total		2093857	100