


HPLC REPORT

Sample: H-10 Analyzed date: 2014-07-16
 Analyst: WGJ
 Lot. No.: P140711-HX594820
 Column: SinoChrom ODS-BP, 4.6*250mm, 5µm
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water

Gradient:

	A	B
0.0min	21%	79%
25.0min	46%	54%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl
 Wavelength: 220nm
 Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	6.524	162960	1.07
2	7.086	229948	1.51
3	7.364	130221	0.8553
4	7.628	599907	3.94
5	7.913	13129418	86.24
6	8.258	355745	2.336
7	8.846	616892	4.052

Total 100