

HPLC REPORT

Sample: K-12

Analyzed date: 2016-01-15

Analyst: ZK

Lot. No.: P160108-HS631192

Column: Kromasil 100-5C18, 4.6*250mm, 5µm

Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile

Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water

Gradient :

	A	B
0.0min	15%	85%
25.0min	40%	60%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl

Wavelength: 220nm

Flow rate: 1.0ml/min

mV
-320

-280

-240

-200

-160

-120

-80

-40

0

-40

-80

2 4 6 8 10 12 14 16 18 20 22 24 min

Rank	Time	Quantity	Area	Height
1	10.193	0.6137	12031	4677
2	10.374	95.9352	1880588	215050
3	10.605	3.4511	67652	5618
Total		100	1960271	225345