


MASS SPECTROMETRY REPORT



Sample Description

Analyzed date: 2016-1-12
Analyst: Xu
Sample: K-15
M.W.: 1765.08
Lot. No.: P160108-SW631325

Instrument

SHIMADZU LCMS-2010EV
Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp.: 250 °C
Block Temp.: 200 °C

Probe Bias: + 4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. Conc.: 50% H₂O/50% ACN