


# HPLC REPORT

Sample: P-16 Analyzed date: 2016-06-11  
 Analyst: LLR  
 Lot. No.: P160604-LL658611  
 Column: Kromasil C18, 4.6\*250mm, 5µm  
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile  
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water  
 Gradient:
 

	A	B
0.0min	12%	88%
25.0min	37%	63%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl  
 Wavelength: 220nm  
 Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	11.767	623603.938	5.7611
2	11.950	9784406.000	90.3920
3	11.950	416406.281	3.8469
Total		10824416.219	100