


MASS SPECTROMETRY REPORT



Sample Description

Analyzed date: 2015-6-11
 Analyst: Han
 Sample: S-12
 M.W.: 1387.53
 Lot. No.: P150604-MX658715

Instrument

SHIMADZU LCMS-2010EV
 Probe: ESI
 Nebulizer Gas Flow: 1.5L/min
 CDL: -20.0v
 CDL Temp.: 250 °C
 Block Temp.: 200 °C

Probe Bias:

+ 4.5kv
 Detector: 1.5kv
 T. Flow: 0.2ml/min
 B. Conc.: 50% H2O/50% ACN