



Sample Description

Analyzed date: 2014/6/15

Analyst: Shen

Sample: I-16

M.W.: 1890.16

Lot. No.: P140605-GY659370

Instrument SHIMADZU LCMS-2020

Probe: ESI Probe Bias: +4.5kv

Nebulizer Gas Flow: 1.5L/min Detector: 1.2kv

CDL: -20.0v T. Flow: 0.2ml/min

CDL Temp.: 250 °C B. Conc.: 50% H<sub>2</sub>O/50% ACN

Block Temp.: 400 °C