

HPLC REPORT

Sample: N-14 Analyzed date: 2014-06-09

Analyst: YHY

Lot. No.: P140601-LL587348

Column: SinoChrom ODS-BP, 4.6*250mm, 5µm

Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile

Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water

Gradient:


	A	B
0.0min	10%	90%
25.0min	35%	65%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl

Wavelength: 220nm

Flow rate: 1.0ml/min

mV



Rank	Time	Area	Conc.
1	8.569	4355	0.1913
2	9.013	2110700	92.69
3	9.328	112976	4.961
4	9.829	49206	2.161
Total		2277237	100