

HPLC REPORT

Sample: T-16

Analyzed date: 2015-08-22

Analyst: ZK

Lot. No.: P150817-SY671689

Column: Kromasil 100-5C18, 4.6*250mm, 5µm

Solvent A A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile

Solvent B B: 0.1% Trifluoroacetic Acid in 100% Water


Gradient :

	A	B
0.0min	10%	90%
25.0min	35%	65%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl

Wavelength: 220nm

Flow rate: 1.0ml/min



Rank	Time	Quantity	Area	Height
1	7.509	6.4140	568190	80106
2	8.168	1.5343	135920	22026
3	8.842	2.0622	182681	43251
4	8.988	89.7813	7953390	1044402
5	9.240	0.2082	18446	9149
Total		100	8858627	1198934