


HPLC REPORT

Sample: S-14 **Analyzed date:** 2014-03-18
Analyst: YHY
Lot. No.: P140228-YS197852
Column: SinoChrom ODS-BP, 4.6*250mm, 5µm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
Gradient:

	A	B
0.0min	9%	91%
25.0min	34%	66%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl
Wavelength: 220nm
Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	13.045	64796.773	1.2536
2	13.610	13687.966	0.2648
3	13.787	4868793.000	94.1931
4	14.048	215015.359	4.1598
5	14.048	6654.337	0.1287
Total		5168947.436	100