


## HPLC REPORT

Sample: N-15  
Analyst: LLR  
Lot. No.: P160604-LL658495  
Column: Kromasil C18, 4.6\*250mm, 5µm  
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile  
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water  
Gradient:

Analyzed date: 2016-06-08

|         | A    | B   |
|---------|------|-----|
| 0.0min  | 14%  | 86% |
| 25.0min | 39%  | 61% |
| 25.1min | 100% | 0%  |
| 30.0min | Stop |     |

Volume: 20µl  
Wavelength: 220nm  
Flow rate: 1.0ml/min



| Rank | Time   | Area    | Conc. |
|------|--------|---------|-------|
| 1    | 11.193 | 98837   | 1.999 |
| 2    | 11.404 | 4751766 | 96.09 |
| 3    | 11.669 | 94551   | 1.912 |

|       |  |         |     |
|-------|--|---------|-----|
| Total |  | 4945154 | 100 |
|-------|--|---------|-----|