




## Sample Description

Analyzed date: 2016-6-8  
Analyst: Han  
Sample: N-15  
M.W.: 1675.86  
Lot. No.: P160604-LL658495

## Instrument

Probe: ESI  
Nebulizer Gas Flow: 1.5L/min  
CDL: -20.0v  
CDL Temp.: 250 °C  
Block Temp.: 200 °C

## SHIMADZU LCMS-2010EV

Probe Bias: + 4.5kv  
Detector: 1.5kv  
T. Flow: 0.2ml/min  
B. Conc.: 50%H<sub>2</sub>O/50%ACN