


HPLC REPORT

Sample: E-15 Analyzed date: 2014-03-10
 Analyst: YHY
 Lot. No.: P140302-YJ566372
 Column: SinoChrom ODS-BP, 4.6*250mm, 5µm
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
 Gradient:

	A	B
0.0min	8%	92%
25.0min	33%	67%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl
 Wavelength: 220nm
 Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	16.663	115041.961	1.1936
2	16.862	9365216.000	97.1693
3	17.380	157785.422	1.6371
Total		9638043.383	100