


# HPLC REPORT

Sample: L-15  
 Analyst: LLR  
 Lot. No.: P140302-YJ566409  
 Column: Kromasil C18, 4.6\*250mm, 5µm  
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile  
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water  
 Gradient:
 

|         | A    | B   |
|---------|------|-----|
| 0.0min  | 18%  | 82% |
| 25.0min | 43%  | 57% |
| 25.1min | 100% | 0%  |
| 30.0min | Stop |     |

Volume: 20µl  
 Wavelength: 220nm  
 Flow rate: 1.0ml/min

Analyzed date: 2014-03-05



| Rank  | Time   | Area        | Conc.   |
|-------|--------|-------------|---------|
| 1     | 15.600 | 2439622.750 | 98.8895 |
| 2     | 15.600 | 27397.047   | 1.1105  |
| Total |        | 2467019.797 | 100     |