


MASS SPECTROMETRY REPORT

15-Jan-2018
11:26:35

Probe: ESI
Cone: 50v
Desolvation Temp :350

Capillary:3.00KV
Extractor: 5v
Gas Flow : 350
Scan ES-
2.36e6

180115-16- 56 (1.035) Cm (56-63)



Sample Description

Analyzed date: 2016-01-15
Analyst: YU
Sample: I-16
M.W.: 1638.79
Lot. No.: P160110-SY630494

Instrument

Waters ZQ2000
Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp.: 250 °C
Block Temp.: 200 °C

Waters ZQ2000

Probe Bias: + 4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. Conc.: 50%H2O/50%ACN