




Sample Description

Analyzed date: 2015-6-20
Analyst: Han
Sample: Q-14
M.W.: 1527.62
Lot. No.: P150604-LL658569

Instrument

SHIMADZU LCMS-2010EV

Probe:

ESI
Nebulizer Gas Flow: 1.5L/min

CDL:

-20.0v

CDL Temp.:

250 °C

Block Temp.:

200 °C

Probe Bias:

+ 4.5kv

Detector:

1.5kv

T. Flow:

0.2ml/min

B. Conc.:

50%H₂O/50%ACN