


HPLC REPORT

Sample: K-14
 Analyst: YHY
 Lot. No.: P140302-YJ358217
 Column: SinoChrom ODS-BP, 4.6*250mm, 5µm
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
 Gradient:

Analyzed date: 2014-03-10

	A	B
0.0min	18%	82%
25.0min	43%	57%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl
 Wavelength: 220nm
 Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	8.790	7089	0.04105
2	8.951	212032	1.228
3	9.548	17049901	98.73
Total		17269022	100