

MASS SPECTROMETRY REPORT


08-Jun-2017
10:26:21

Probe: ESI
Cone: 50v
Desolvation Temp :350

Capillary:3.00KV
Extractor: 5v
Gas Flow : 350

Scan ES+
1.90e7

170608-13 462 (8.539)



Sample Description

Analyzed date: 2014-06-08
Analyst: YU
Sample: R-13
M.W.: 1523.74
Lot. No.: P140602-CJ587837

Instrument

Waters ZQ2000
Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp.: 250 °C
Block Temp.: 200 °C

Waters ZQ2000

Probe Bias: + 4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. Conc.: 50%H2O/50%ACN