

MASS SPECTROMETRY REPORT


11-Jan-2018
10:31:29

Probe: ESI
Cone: 50v
Desovation Temp :350

Capillary:3.00KV
Extractor: 5v
Gas Flow : 350

Scan ES+
3.92e6

180111-20 182 (3.364) Cm (179:189)



Sample Description

Analyzed date: 2016-01-11
Analyst: YU
Sample: N-16
M.W.: 1944.12
Lot. No.: P160109-SY630463

Instrument

Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp.: 250 °C
Block Temp.: 200 °C

Waters ZQ2000

Probe Bias: + 4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. Conc.: 50%H2O/50%ACN