

MASS SPECTROMETRY REPORT


19-Jun-2017
14:40:10

Probe: ESI
Cone: 50v
Desovation Temp :350

Capillary:3.00KV
Extractor: 5v
Gas Flow : 350

Scan ES-
8.60e3

170619-29- 123 (2.273) Cm (113:127)



Sample Description

Analyzed date: 2014-06-19
Analyst: YU
Sample: A-14
M.W.: 1468.54
Lot. No.: P140615-SY587874

Instrument

Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
CDL Temp.: 250 °C
Block Temp.: 200 °C

Waters ZQ2000

Probe Bias: + 4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. Conc.: 50%H2O/50%ACN