


HPLC REPORT

Sample: C-16
 Analyst: WGJ
 Lot. No.: P140705-YJ593517
 Column: SinoChrom ODS-BP, 4.6*250mm, 5µm
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
 Gradient:

Analyzed date: 2014-07-12

	A	B
0.0min	24%	76%
25.0min	49%	51%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl
 Wavelength: 220nm
 Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	7.401	147860	1.231
2	7.919	687538	5.723
3	8.180	362811	3.02
4	8.623	10264534	85.44
5	8.870	551025	4.587

Total 100