


HPLC REPORT

Sample: E-16 Analyzed date: 2016-04-22
Analyst: WGJ
Lot. No.: P160411-GY649104
Column: SinoChrom ODS-BP, 4.6*250mm, 5µm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water

Gradient:

	A	B
0.0min	19%	81%
25.0min	44%	56%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl
Wavelength: 220nm
Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	9.510	40611	1.065
2	9.619	3497905	91.72
3	10.084	275218	7.216
Total			100