


MASS SPECTROMETRY REPORT

Inten. (x 100,000)



Sample Description

Analyzed date: 2016-4-24
 Analyst: Han
 Sample: E-16
 M.W.: 1652.81
 Lot. No.: P160411-GY649104

Instrument

Probe: ESI
 Nebulizer Gas Flow: 1.5L/min
 CDL: -20.0v
 CDL Temp.: 250 °C
 Block Temp.: 200 °C

SHIMADZU LCMS-2010EV

Probe Bias: + 4.5kv
 Detector: 1.5kv
 T. Flow: 0.2ml/min
 B. Conc.: 50% H2O/50% ACN