


HPLC REPORT

Sample: T-14
Analyst: LLR
Lot. No.: P140705-YS593580
Column: Kromasil C18, 4.6*250mm, 5µm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
Gradient:
0.0min 14% 86%
25.0min 39% 61%
25.1min 100% 0%
30.0min Stop
Volume: 20µl
Wavelength: 220nm
Flow rate: 1.0ml/min

Analyzed date: 2014-07-16



Rank	Time	Area	Conc.
1	10.148	13053.017	0.3728
2	12.123	46942.098	1.3405
3	12.272	3107874.500	88.7530
4	12.272	296736.125	8.4740
5	27.710	37106.730	1.0597
Total		3501712.470	100