


HPLC REPORT

Sample: K-11 Analyzed date: 2016-04-17
Analyst: WGJ
Lot. No.: P160412-HX649172
Column: SinoChrom ODS-BP, 4.6*250mm, 5µm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
Gradient:

	A	B
0.0min	4%	96%
25.0min	29%	71%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl
Wavelength: 220nm
Flow rate: 1.0ml/min
mV



Rank	Time	Area	Conc.
1	8.716	30758	0.7832
2	9.462	30705	0.7818
3	9.662	3568700	90.86
4	9.992	297342	7.571
Total			100