


HPLC REPORT

Sample: Y-16
 Analyst: LLR
 Lot. No.: P140705-YS593592
 Column: Kromasil C18, 4.6*250mm, 5µm
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
 Gradient:

Analyzed date: 2014-07-10

	A	B
0.0min	18%	82%
25.0min	43%	57%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl
 Wavelength: 220nm
 Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	8.923	19961	0.8299
2	9.251	47174	1.961
3	9.508	42063	1.749
4	9.655	2058894	85.6
5	9.847	163953	6.816
6	10.216	73288	3.047
Total		2405333	100