


HPLC REPORT

Sample: Y-14
Analyst: WGJ
Lot. No.: P140705-YS593620
Column: SinoChrom ODS-BP, 4.6*250mm, 5µm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
Gradient:
0.0min 18% 82%
25.0min 43% 57%
25.1min 100% 0%
30.0min Stop
Volume: 20µl
Wavelength: 220nm
Flow rate: 1.0ml/min

Analyzed date: 2014-07-13



Rank	Time	Area	Conc.
1	11.347	25698.750	12.2883
2	11.723	1804178.125	86.2765
3	11.723	30010.721	1.4351
Total			100