

HPLC REPORT

Sample: N-16-2
 Analyst: WGJ
 Lot. No.: P160509-HX654162
 Column: SinoChrom ODS-BP, 4.6*250mm, 5µm
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water


Analyzed date: 2016-05-15

Gradient:

	A	B
0.0min	11%	89%
25.0min	36%	64%
25.1min	100%	0%
30.0min		

Stop

Volume: 20µl
 Wavelength: 220nm
 Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	11.985	14902	0.4876
2	12.128	2637501	86.29
3	12.435	403991	13.22

Total 100