


HPLC REPORT

Sample: D-12
Analyst: HXH
Lot. No.: P140714-SY546400
Column: Gemini-NX 5 μ C18 110A, 4.6*250mm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
Gradient:

Analyzed date: 2014-07-17

| | A | B |
|---------|------|-----|
| 0.0min | 15% | 85% |
| 25.0min | 40% | 60% |
| 25.1min | 100% | 0% |
| 30.0min | Stop | |

Volume: 20 μ l
Wavelength: 220nm
Flow rate: 1.0ml/min



| Rank | Time | Conc. | Area | Height |
|-------|--------|---------|----------|--------|
| 1 | 11.981 | 3.5803 | 535233 | 18094 |
| 2 | 12.393 | 93.8929 | 14036206 | 726839 |
| 3 | 12.788 | 0.7995 | 119519 | 18005 |
| 4 | 13.678 | 1.7273 | 258215 | 14383 |
| Total | 100 | | 14949173 | 777321 |