


MASS SPECTROMETRY REPORT



Sample Description

Analyzed date: 2014-6-9
 Analyst: Long
 Sample: V-14
 M.W.: 1582.82
 Lot. No.: P140601-LL294919

Instrument

SHIMADZU LCMS-2010EV
 Probe: ESI
 Nebulizer Gas Flow: 1.5L/min
 CDL: -20.0v
 CDL Temp.: 250 °C
 Block Temp.: 200 °C

Probe Bias: + 4.5kv
 Detector: 1.5kv
 T. Flow: 0.2ml/min
 B. Conc.: 50% H2O/50% ACN