


HPLC REPORT

Sample: K-13
 Analyst: YW
 Lot. No.: P160110-SY630457
 Column: 4.6×250mm, Kromasil 100-5-C18
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
 Gradient:

	A	B
0.0min	12%	88%
25.0min	37%	63%
25.1min	100%	0%
30.0min	Stop	

Volume: 5µl
 Wavelength: 220nm
 Flow rate: 1.0ml/min

Analyzed date: 2016-01-12



Peak	Time	Height	Area	Conc.
1	12.348	5817.985	39592.215	0.7949
2	12.508	3258.711	21872.746	0.4391
3	12.727	7706.092	48659.266	0.9769
4	12.878	533007.625	4557175.000	91.4901
5	12.878	2857.033	16049.507	0.3222
6	13.540	30342.752	205064.859	4.1169
7	13.843	11280.813	74727.313	1.5002
8	15.998	2814.637	17916.502	0.3597
Total				100.0000