


HPLC REPORT

Sample: S-13 Analyzed date: 2014-06-07
Analyst: LLR
Lot. No.: P140601-LL118022
Column: Kromasil C18, 4.6*250mm, 5µm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
Gradient:
0.0min 22% 78%
25.0min 47% 53%
25.1min 100% 0%
30.0min Stop
Volume: 20µl
Wavelength: 220nm
Flow rate: 1.0ml/min



| Rank | Time | Area | Conc. |
|-------|--------|-------------|---------|
| 1 | 11.035 | 84581.047 | 4.2770 |
| 2 | 11.213 | 1760921.375 | 89.0451 |
| 3 | 11.213 | 22321.443 | 1.1287 |
| 4 | 12.933 | 109736.953 | 5.5491 |
| Total | | 1977560.818 | 100 |