


MASS SPECTROMETRY REPORT

Inten. (x 10,000,000)



Sample Description

Analyzed date: 2014-6-6
 Analyst: Long
 Sample: S-13
 M.W.: 1356.47
 Lot. No.: P140601-LL118022

Instrument

Probe: ESI
 Nebulizer Gas Flow: 1.5L/min
 CDL: -20.0v
 CDL Temp.: 250 °C
 Block Temp.: 200 °C

SHIMADZU LCMS-2010EV

Probe Bias: + 4.5kv
 Detector: 1.5kv
 T. Flow: 0.2ml/min
 B. Conc.: 50% H2O/50% ACN