

# HPLC REPORT

Sample: P-16 Analyzed date: 2015-05-24

Analyst: RZG

Lot. No.: P150511-GY654305

Column: SinoChrom ODS-BP, 4.6\*250mm, 5µm

Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile


Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water

Gradient:	A	B
0.0min	10%	90%
25.0min	35%	65%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl

Wavelength: 220nm

Flow rate: 1.0ml/min



Rank	Time	Conc	Area
1	16.300	0.7872	64492
2	16.744	93.96	7698049
3	17.608	5.252	430253
Total		100	8192794