


MASS SPECTROMETRY REPORT



Sample Description

Analyzed date: 2016-5-18
 Analyst: SHAO
 Sample: G-16
 M.W.: 1530.95
 Lot. No.: P160508-MJ654319

Instrument

SHIMADZU LCMS-2020
 Probe: ESI
 Nebulizer Gas Flow: 1.5L/min
 CDL: -20.0v
 CDL Temp.: 250 °C
 Block Temp.: 200 °C

SHIMADZU LCMS-2020

Probe Bias: + 4.5kv
 Detector: 1.5kv
 T. Flow: 0.2ml/min
 B. Conc.: 50% H2O/50% ACN