


HPLC REPORT

Sample: L-16
Analyst: WJ
Lot. No.: P140724-HS598231
Column: Gemini-NX 5 μ C18 110A, 4.6*250mm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
Gradient:

Analyzed date: 2014-8-8

	A	B
0.0min	20%	80%
25.0min	80%	20%
25.1min	100%	0%
30.0min	Stop	

Volume: 20 μ l
Wavelength: 220nm
Flow rate: 1.0ml/min



Rank	Time	Conc.	Area	Height
1	10.304	3.0435	261909	56056
2	10.447	96.0675	8267005	907918
3	10.704	0.8890	76504	22737
Total		100	8605418	986711