

HPLC REPORT

Sample: D-16 Analyzed date: 2016-04-21

Analyst: WGJ

Lot. No.: P160411-GY649133

Column: SinoChrom ODS-BP, 4.6*250mm, 5µm

Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile

Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water

Gradient:


	A	B
0.0min	8%	92%
25.0min	33%	67%
25.1min	100%	0%
30.0min	Stop	

Volume: 20µl

Wavelength: 220nm

Flow rate: 1.0ml/min

mV



Rank	Time	Area	Conc.
1	8.604	95988	2.998
2	8.730	2909644	90.87
3	8.922	196367	6.133
Total			100